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Segreteria Tel. 0471-438306	Sekretariat Tel. 0471-438306
Riferimenti bibliografici [RIF.1] WT1 expression as a marker of minimal residual disease predicts outcome in acute myeloid leukemia when measured post-consolidation. Leukemia Research, 2012,36:453-458 [RIF.2] Real-time assessment of relapse risk based on the WT1 marker in acute leukemia and myelodysplastic syndrome patients after hematopoietic cell transplantation. BMT, 2015, 50: 26-33 [RIF.3] Separatori di linfociti Arnika: foglio illustrativo [RIF.4] Lymphoprep: foglio illustrativo [RIF.5] Maxwell CSC RNA Blood Kit TM434 Italian v.10_22 [RIF.6] Manuale del kit ipsogen® RT [RIF.7] Manuale del kit ipsogen® WT1 [RIF.8] Role of Sequential Monitoring of WT1 Expression as a Reliable Marker for Prognosis and Minimal Residual Disease Monitoring in Acute Myeloid Leukemia With Normal Karyotype Patients. Clinical Lymphoma, Myeloma & Leukemia, 2017,17: e521-527 [RIF.9] WT1 gene expression: an excellent tool for monitoring minimal residual disease in 70% of acute myeloid leukaemia patients –results from a single-centre. Bjh, 2004, 125:590-600 [RIF.10] Leukaemia relapse after allogeneic transplants for acute myeloid leukemia: predictive role of WT1 expression. Bjh, 2013,160: 503-509 [RIF.11] Role of Sequential Monitoring of WT1 Gene Expression in Patients with Acute Myeloid Leukemia for the Early Detection of Leukemia Relapse. Clinical Lymphoma, Myeloma & Leukemia, 2017,17: e521-527 [RIF. 12] Diagnosis and management of AML in adults: 2022 recommendations from an international expert panel on behalf of the ELN, Blood, 2022, 140: 1345-1377	Literatur [RIF.1] WT1 expression as a marker of minimal residual disease predicts outcome in acute myeloid leukemia when measured post-consolidation. Leukemia Research, 2012,36:453-458 [RIF.2] Real-time assessment of relapse risk based on the WT1 marker in acute leukemia and myelodysplastic syndrome patients after hematopoietic cell transplantation. BMT, 2015, 50: 26-33 [RIF.3] Separatori di linfociti Arnika: foglio illustrativo [RIF.4] Lymphoprep: foglio illustrativo [RIF.5] Maxwell CSC RNA Blood Kit TM434 Italian v.10_22 [RIF.6] Manuale del kit ipsogen® RT [RIF.7] Manuale del kit ipsogen® WT1 [RIF.8] Role of Sequential Monitoring of WT1 Expression as a Reliable Marker for Prognosis and Minimal Residual Disease Monitoring in Acute Myeloid Leukemia With Normal Karyotype Patients. Clinical Lymphoma, Myeloma & Leukemia, 2017,17: e521-527 [RIF.9] WT1 gene expression: an excellent tool for monitoring minimal residual disease in 70% of acute myeloid leukaemia patients –results from a single-centre. Bjh, 2004, 125:590-600 [RIF.10] Leukaemia relapse after allogeneic transplants for acute myeloid leukemia: predictive role of WT1 expression. Bjh, 2013,160: 503-509 [RIF.11] Role of Sequential Monitoring of WT1 Gene Expression in Patients with Acute Myeloid Leukemia for the Early Detection of Leukemia Relapse. Clinical Lymphoma, Myeloma & Leukemia, 2017,17: e521-527 [RIF. 12] Diagnosis and management of AML in adults: 2022 recommendations from an international expert panel on behalf of the ELN, Blood, 2022, 140: 1345-1377
Aggiornato il 27/03/2026	Aktualisiert am 27/03/2026
Prossimo aggiornamento 27/03/2029	Nächste Aktualisierung am 27/03/2029