

HemaTrack-ALL: Taking advantage of Unique Molecular Identifiers and NGS for accurately monitoring MRD in ALL


U.-P. Guenther¹, V. Surendranath¹, C. Pahlke¹, M. Pahlke¹, T. Lingl¹, T. Kasper¹, P. Dasgupta², C. Eckert³, V. Saha², V. Lange¹

¹ DKMS Life Science Lab gGmbH, Dresden, Germany
² TATA Translational Cancer Research Centre, Kolkata, India
³ Charité – Universitaetsmedizin, Berlin, Germany

Introduction

Measurable Residual Disease (MRD) is the single best predictor for survival outcomes in Acute Lymphoblastic Leukemia (ALL). PCR-based MRD assays have emerged as gold standard in Europe promising higher specificity and reproducibility than flow cytometry-based approaches. The use of NGS may outperform laborious patient-specific assay development in qPCR-based quantification approaches, as it provides high level of specificity due to single nucleotide resolution while delivering sufficient data to achieve high sensitivity. However, biases in PCR amplification may severely distort accuracy of quantification. To address this issue, we developed a robust NGS-based approach for MRD monitoring in ALL, that overcomes existing limitations by incorporating Unique Molecular Identifiers (UMIs).

Methods




Schematic workflow of HemaTrack-ALL. HemaTrack amplifies clonotype sequences in three major steps whereas each step is followed by a clean up procedure:

- 1. UMI-tagging step:** UMIs are appended with two cycles of PCR containing 500ng of genomic DNA in two highly multiplexed reactions targeting IGH, IGK, TRG and TRD with primer set 1 and optionally, TRB with primer set 2. A third reaction amplifies a reference gene in a singleplex reaction. Noteworthy, the same workflow is applied to patient samples at diagnosis or at a follow-up timepoint except that reactions containing primer set 1 or 2 are performed in triplicates in case of follow-up samples.
- 2. Index PCR:** Universal primers containing sample barcodes and Illumina adapter sequences are added and a second PCR is performed for a limited number of cycles. Triplicate reaction are pooled at this step in case of samples from follow-up timepoints
- 3. Final PCR:** Primers containing only Illumina adapters are added to the purified products of the Index PCR.


Results

Assay specifications

Highly multiplexed PCR workflow




Accurate quantification




HemaTrack-ALL quantifies accurately. HemaTrack-ALL was performed on genomic DNA of several ALL cell lines with known clonotypes and frequencies were diluted serially into blood DNA from healthy individuals.

High sensitivity




Comparison with gold standard

Patient cohort and samples



HemaTrack-ALL monitors more clonotypes



Good concordance with RQ-PCR in MRD calls


		RQ-PCR	
		MRD	
HemaTrack-ALL	+	7	5
	-	1	15
		88%	75%
		sensitivity	specificity

HemaTrack-ALL achieves high level of sensitivity and specificity in comparison to gold standard (RQ-PCR). 2x2 contingency table comparing the MRD statuses from 28 follow-up timepoints called by both methods.

Conclusion

HemaTrack-ALL is an innovative NGS-based, end-to-end workflow which features the identification of markers of ALL blasts and their subsequent high sensitivity quantitative tracking. The assay presents with a straightforward way of quantification, meets the clinical requirements in terms of accuracy and sensitivity and will be further evaluated in real-live clinical routine in a prospective clinical trial in India next year.

Outlook



Acknowledgements

Team TTCRC

Team Charité

Team LSL

Anja Klußmeier	Lisa-Marie Moskwa	Patrick Paul
Daniela Krater	Henrike Indrischek	Rajath TS
Daniela Walter	Bianca Schöne	
Doreen Tolksdorf	Gerhard Schöfl	
Geoffrey Behrens		
Thomas Schäfer		

