

phits

PRODUCT PORTFOLIO

FEATURED DESIGN PARAMETERS

	EVERYDAY USE			SPORTS				WORKWEAR		
	NARROW	COMFORT	WIDE	RUNNING	SOCCER/ FOOTBALL	GOLF	CYCLING	ALPINE SKI	NORDIC SKI	SAFETY
										
Dynamic footscan analysis protocol	walking			walking + running		walking	walking + running + one legged forefoot squat	walking+ running one legged forefoot squat or double poling simulation		walking
Standard top cover	3 mm EVA (30, 35, 40 shore)			6 mm EVA (30,35, 40 shore)		3 mm EVA+ synth leather (30, 35, 40 shore)	2 mm PU soft	2 mm PU soft		3 mm EVA carbon (40 shore, antistatic combination)
Other top cover options	6 mm EVA (30, 35, 40 shore) 3 mm PU soft (20 shore) 6 mm PU soft (20 shore)			3 mm EVA (30, 35, 40 shore) 3 mm PU soft (20 shore) 6 mm PU soft (20 shore)		6 mm EVA (30, 35, 40 shore) 3 mm PU soft (20 shore) 6 mm PU soft (20 shore)	2 mm EVA	2 mm EVA		6 mm EVA carbon
Synth leather top cover available	yes			yes		yes	yes	yes		yes
'Normal-slim' and 'Ortho-slim' Phits option (no lattice structure)	yes			yes		yes	no	no		yes
Ortho type available	yes			yes		yes	yes	yes		yes
Directional stiffness prediction (adjustable -30°/+30°)	yes*			yes*		yes*	yes	yes		yes*
Local stiffness differentiation (adjustable)	yes			yes		yes	yes	yes		yes
Additional heel cushioning	plantar fasciitis, heel spur, fat pad			plantar fasciitis, heel spur, fat pad		plantar fasciitis, heel spur, fat pad	no	no		plantar fasciitis, heel spur, fat pad
Low heel cup option	yes			yes		yes	yes	yes		yes
Heel wedge	range from -5° up to 7°*			range from -5° up to 7° *		range from -5° up to 7° *	range from -5° up to 7°	range from -5° up to 7°		range from -5° up to 7°*
Heel offset	1-6 mm*			1-6 mm*		1-6 mm*	1-6mm	1-6mm		1-6 mm*
Low lateral edge option	yes			yes		yes	yes	yes		yes
Meta bar available	1-3 mm			1-3 mm		1-3 mm	no	1-3 mm		1-3 mm
Meta pads available	teardrop/T-form 1-5 mm			teardrop/T-form 1-5 mm		teardrop/T-form 1-5 mm	teardrop/T-form 1-5 mm	teardrop/T-form 1-5 mm		teardrop/T-form 1-5 mm
Forefoot correction	Meta 1, Meta 1-2, Meta 2-4, Meta 2-5, Meta 5 height 2-4 mm Flexible/rigid*			Meta 1, Meta 1-2, Meta 2-4, Meta 2-5, Meta 5 height 2-4 mm Flexible/rigid*		Meta 1, Meta 1-2, Meta 2-4, Meta 2-5, Meta 5 height 2-4 mm Flexible/rigid*	Complete forefoot medial/lateral 1-4 mm rigid	Complete forefoot medial/lateral 1-4 mm Flexible		Meta 1, Meta 1-2, Meta 2-4, Meta 2-5, Meta 5 height 2-4 mm Flexible/rigid*
Patient name & serial number	yes			yes		yes	yes	yes		yes
* not available for Slim Phits										